

ABI-600CL

Automated Single-level, Multi-level, & Stress Testing
CPT 93922, 93923, 93924

Newman Medical

simpleABI Systems from Newman Medical use the power of computers to produce the reimbursement reports you need quickly, easily, and accurately. The easy to use PC-based systems' minimal learning curve makes it easy to begin performing exams. Automatic index calculations make the simpleABI system make exams as easy as taking blood pressures. **Print to any standard office printer or save and load directly to EMR.** The diagnosis is made using the gold standard Doppler and simple pulse volume recording (PVR) waveforms. The ABI-600CL system meets all requirements, including dual ankle, for CPT code 93922, 93923, & 93924.

ABI-600CL Features:

- **Cuff-Link™ automated cuff select, inflate, calculate**
- PC-based ABI system
 - Easy staff training - simple to learn and use
 - Point & Click reporting right on the computer
- Uses "Gold Standard" 8MHz Doppler
- Easy "Click & Capture" PVR waveforms
- Customize reports with practice information and logo
- **Print to any standard office printer**
- **Attach multiple formats directly to EMR records**
- Full 2-year product warranty

ABI-600CL Kit Contains:

- Cuff-Link™-8
- Single Level ABI & TBI Testing
- Multi-level Segmental Testing
- **ABI Stress Testing**
- Sensitive 8MHz Doppler
- PPG Probe
- 9 cuffs (4x10cm, 4x12cm, 1x1.9cm)



ABI-600CL Options:

- Computer
- Roll Stand
- Printer
- Mobile Carry Bag

Distributed By:



At the Heart of Innovative Technology Solutions

Contact us today to place an order

Phone: 888-LynnMed / 888-596-6633

Fax (248) 560-0500

www.lynnmed.com



NM0004B